


























## Applicable symbology according to ISO 15223-1

Category	Information to supply	Graphic	
Generic product identification	Medical device name	n.a	
	Registered trademark		
	Indication of medical device		
	Identification of the package contents	Composition	Cotton 
			Polyester/Viscose 
		Dimension 	
		Plies 	
		Gauze threads/cm <sup>2</sup> 	
		Non-Woven grammage 	
		Number of units 	
	Product identification	Catalogue number	
	Production data	Unique device identifier	
		Batch Code	
		SN Number	
		Use-by-date	
	Country of origin		
Manufacturer	Name		
	Address		
	Phone number		
	Web address		
CE marking	CE marking with notified body number		
Indication for single use	Do not re-use		
Sterilization method	Sterilized using steam or EtO		
Indication for proper use	Do not resterilize		
Special handling conditions	Keep away from sunlight		
	Keep dry		
	Do not use if package is damaged and consult instructions for use		
Additional info	Recyclable packaging		
Website	Company website		
Reference to the absence of latex	Symbol "Latex free"	LATEX FREE   SIN LATEX   SANS LATEX   SEM LATEX	
Warnings	Warning to verify the number of pieces		
	Warning to use in its final fold	n.a	
Distributor (when applicable)	Name		
	Address		
	Phone number	